



MEDICAL PROFESSIONALS'
TRUSTED BRAND

UROLOGY

PRODUCT LIST

GUIDEWIRES (5-12)

- PTFE Guidewire
- Twister Guidewire
- Hydro Twister™ Guidewire
- Tracer™ Guidewire
- Amplatz Extra Stiff Guidewire
- Amplatz Ultra Stiff™ Guidewire
- Prime Wire™

STENTS (13-23)

- P Flex Stent With String (Ureteral Stent / Double J Stent)
- P Flex Stent Set With String (Ureteral Stent Set / Double J Stent Set)
- Endopyelotomy Stent
- Urethral Stent

URETERAL ACCESS SHEATH (24-25)

- Praveza® Ureteral Access Sheath

CATHETERS (26-35)

- Ureteral Catheter
- Mono J Catheter
- PCN Catheter
- PCN Catheter with Locking Mechanism
- PCN Catheter with Needle
(Nephrostomy Drainage Set Direct Puncture)
- Malecot Catheter
- Re-entry Malecot Drainage Catheter
- Suprapubic Catheter (SPC)
- CI Catheter

DRAINAGE SETS (36-40)

- PCN Set - Pigtail
- PCN Set - Pigtail with Locking Mechanism
- PCN Set - Malecot
- Suprapubic Drainage Catheter Set (SDS)

STONE BASKETS (41-48)

- N Quadra™ - Nitinol Stone Basket With Tip
- N Quadra™ Z - Nitinol Stone Basket Zero Tip
- N Quadra™ Helical - Nitinol Stone Basket Helical With Tip
- N Quadra™ Helical Z - Nitinol Stone Basket Helical Zero Tip
- N Quadra™ Flexible - Nitinol Stone Basket Flexible
- Delta Basket
- PCNL Basket

PRODUCT LIST

DILATORS (49-63)

- Fascial Dilator Set
- Ureteral Dilator Set
- Single Step™ Ureteral Dilator
- Urethral Dilator Set
- Urethral Dilator Set - S Curved
- Urethral Self Dilator
- Amplatz Sheath
- Amplatz Renal Dilator
- Amplatz Renal Dilation Set - S
- Amplatz Renal Dilation Set - F
- Mini Perc Renal Dilation Set
- Ureteral Balloon Dilator
- Nephrostomy Balloon Dilator With Amplatz Sheath
- Screw Dilator

NEEDLES (64-66)

- Initial Puncture Needle - Two Part
- Sure Puncture™ Needle (Initial Puncture Needle - Three Part)

BIOPSY (67-68)

- Biopsy Gun - Semi Automatic

ACCESSORIES (69-75)

- Meatal Dilator - Male
- Dalela Female Self Urethral Dilator
- Bladder Evacuator
- Penile Clamp
- Irrigation Set (TUR Set)
- Cath Tap

GUIDEWIRES



- PTFE Guidewire
- Twister Guidewire
- Hydro Twister™ Guidewire
- Tracer™ Guidewire
- Amplatz Extra Stiff Guidewire
- Amplatz Ultra Stiff™ Guidewire
- Prime Wire™

UROLOGY

PTFE GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures (URS).

Key Features

- PTFE Coating for smooth insertion
- Atraumatic floppy tip for less trauma
- Uniform outer diameter

Types


- J Tip
- Straight Tip

Size (in inch)

0.018, 0.025, 0.028, 0.032, 0.035
6 0.038

Length (in cm)

150



PTFE GUIDEWIRE - J TIP

Product Code	Size (in inch)	Length (in cm)
BPTSO18I50BJT	0.018	150
BPTSO25I50BJT	0.025	150
BPTSO28I50BJT	0.028	150
BPTSO32I50BJT	0.032	150
BPTSO35I50BJT	0.035	150
BPTSO38I50BJT	0.038	150

PTFE GUIDEWIRE - STRAIGHT TIP

Product Code	Size (in inch)	Length (in cm)
BPTSO18I50BST	0.018	150
BPTSO25I50BST	0.025	150
BPTSO28I50BST	0.028	150
BPTSO32I50BST	0.032	150
BPTSO35I50BST	0.035	150
BPTSO38I50BST	0.038	150

TWISTER GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices (URS procedure, for ureteral stenting).

Key Features

- Is constructed with kink resistant nitinol core with appropriate stiffness for deployment of devices such as stents during a procedure
- Atraumatic floppy tip of 5 cm facilitates easy insertion and access
- Blue & White sleeve provides excellent visualization and easy maneuvering during diagnostic and interventional procedure

Type


Straight Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150



TWISTER GUIDEWIRE - STRAIGHT TIP

Product Code	Size (in inch)	Length (in cm)
BTWN025I50BST	0.025	150
BTWN032I50BST	0.032	150
BTWN035I50BST	0.035	150
BTWN038I50BST	0.038	150

HYDRO TWISTER™ GUIDEWIRE

Used in placement of medical devices while gaining easy access to desired anatomical location (URS & RIRS).

Key Features

- Hydrophilic coated tip and blue & white sleeve provides easy, low friction advancement with excellent visualization and maneuverability during diagnostic and interventional procedures
- Atraumatic floppy tip of 5 cm facilitates easy insertion and access
- Is constructed with kink resistant nitinol core with appropriate stiffness for deployment of devices such as stents during a procedure

Advantages

- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy
- Serves two purpose - access and stenting with single wire

Type


Straight Tip

Size (in inch)

0.025, 0.032, 0.035 & 0.038

Length (in cm)

150



HYDRO TWISTER™ GUIDEWIRE - STRAIGHT TIP

Product Code	Size (in inch)	Length (in cm)
BHTN025150BST	0.025	150
BHTN032150BST	0.032	150
BHTN035150BST	0.035	150
BHTN038150BST	0.038	150

TRACER™ GUIDEWIRE (HYDROPHILIC GUIDEWIRE)

Used to gain smooth and easy access to desired anatomical location and to assist in placement of medical devices.

Key Features

- Lubricious hydrophilic coating provides easy passage in tortuous anatomy and through lesions
- Radiopaque polymer jacket enhances fluoroscopic visualization
- Kink resistant nitinol core for optimal performance

Types

- Straight Tip
- Angled Tip

Size (in inch)

0.018, 0.025, 0.032, 0.035 & 0.038

Length (in cm)

150



TRACER™ GUIDEWIRE - ANGLED TIP

Product Code	Size (in inch)	Length (in cm)
BHWN018150BAT	0.018	150
BHWN025150BAT	0.025	150
BHWN032150BAT	0.032	150
BHWN035150BAT	0.035	150
BHWN038150BAT	0.038	150

TRACER™ GUIDEWIRE - STRAIGHT TIP

Product Code	Size (in inch)	Length (in cm)
BHWN018150BST	0.018	150
BHWN025150BST	0.025	150
BHWN032150BST	0.032	150
BHWN035150BST	0.035	150
BHWN038150BST	0.038	150

AMPLATZ EXTRA STIFF GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices and heavier instruments (PCN & PCNL).

Key Features

- PTFE coating for smooth insertion
- Atraumatic floppy tip for less trauma
- The stiffness of the wire provides support for the advancement of devices such as dilators, catheters and instruments

Type

- Straight Tip
- J Tip

Size (in inch)

0.032 & 0.035

Length (in cm)

80



AMPLATZ EXTRA STIFF GUIDEWIRE - J TIP

Product Code	Size (in inch)	Length (in cm)
BAES032080BJT	0.032	80
BAES035080BJT	0.035	80

AMPLATZ EXTRA STIFF GUIDEWIRE - STRAIGHT TIP

Product Code	Size (in inch)	Length (in cm)
BAES032080BST	0.032	80
BAES035080BST	0.035	80

AMPLATZ ULTRA STIFF™ GUIDEWIRE

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures (PCNL).

Key Features

- PTFE coated stainless steel mandril with 5 cm atraumatic floppy tip
- Smooth transition from stiff mandril to the flexible tip
- The stiffness of the mandril design provides support for the advancement of devices such as dilators, catheters and instruments

Type

J Tip

Size (in inch)

0.035

Length (in cm)

80



AMPLATZ ULTRA STIFF™ GUIDEWIRE - J TIP

Product Code	Size (in inch)	Length (in cm)
BASSO35080BJT	0.035	80

PRIME WIRE™

NEW LAUNCH (Gold Standard Guidewire for Urology)

Used to gain access to desired anatomical location and to assist in placement of medical devices during diagnostic or interventional procedures.

Key Features

- Designed to combine the benefits of multiple specialty guidewires in one guidewire for urological procedures
- 5 cm hydrophilic floppy tip provides easy passage in tortuous anatomy
- Teflon coating provides for low friction navigation and easy advancement of devices

Advantages

- One wire for all urology procedures
- Serves multiple purposes - access and stenting with single wire
- Avoid exchange of wires
- Less procedural time
- Easy access in tortuous anatomy

Types


- Straight Tip
- Angled Tip

Size (in inch)

0.032, 0.035 & 0.038

Length (in cm)

150



PRIME WIRE™ - STRAIGHT TIP			PRIME WIRE™ - ANGLED TIP		
Product Code	Size (in inch)	Length (in cm)	Product Code	Size (in inch)	Length (in cm)
BPTN032150BHS	0.032	150	BPTN032150BHA	0.032	150
BPTN035150BHS	0.035	150	BPTN035150BHA	0.035	150
BPTN038150BHS	0.038	150	BPTN038150BHA	0.038	150

STENTS

- P Flex Stent With String (Ureteral Stent / Double J Stent)
- P Flex Stent Set With String (Ureteral Stent Set / Double J Stent Set)
- Endopyelotomy Stent
- Urethral Stent

UROLOGY

P FLEX STENT WITH STRING

(URETERAL STENT / DOUBLE J STENT)

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

Key Features


- For indwelling up to 30 days
- Polyurethane based stent has uniform diameter and smooth surface reducing encrustation
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization
- Stent with string which helps in repositioning

Set Contains

Stent with string, pusher & clamp

Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- With or Without Hydrophilic Coating



Standard

Size (in Fr)	Length (in cm)
4.7, 4.8, 5, 6, 7 & 8	12, 14, 16, 18, 20, 22, 24, 26 28 & 30

Multiloop

Size (in Fr)	Length (in cm)
4.7, 4.8, 5, 6, 7 & 8	18-24 & 24-30

Pediatric

Size (in Fr)	Length (in cm)
3 & 4	10, 12, 14 & 16

P FLEX STENT BOTH ENDS OPEN
WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B1	3	10
BUSB003012B1	3	12
BUSB003014B1	3	14
BUSB003016B1	3	16
BUSB004010B1	4	10
BUSB004012B1	4	12
BUSB004014B1	4	14
BUSB004016B1	4	16
BUSB047012B1	4.7	12
BUSB047014B1	4.7	14
BUSB047016B1	4.7	16
BUSB047018B1	4.7	18
BUSB047020B1	4.7	20
BUSB047022B1	4.7	22
BUSB047024B1	4.7	24
BUSB047026B1	4.7	26
BUSB047028B1	4.7	28
BUSB047030B1	4.7	30
BUSB048012B1	4.8	12
BUSB048014B1	4.8	14
BUSB048016B1	4.8	16
BUSB048018B1	4.8	18
BUSB048020B1	4.8	20
BUSB048022B1	4.8	22
BUSB048024B1	4.8	24
BUSB048026B1	4.8	26
BUSB048028B1	4.8	28
BUSB048030B1	4.8	30
BUSB005012B1	5	12
BUSB005014B1	5	14
BUSB005016B1	5	16
BUSB005018B1	5	18
BUSB005020B1	5	20
BUSB005022B1	5	22
BUSB005024B1	5	24
BUSB005026B1	5	26
BUSB005028B1	5	28
BUSB005030B1	5	30
BUSB006012B1	6	12
BUSB006014B1	6	14
BUSB006016B1	6	16
BUSB006018B1	6	18
BUSB006020B1	6	20
BUSB006022B1	6	22
BUSB006024B1	6	24
BUSB006026B1	6	26
BUSB006028B1	6	28
BUSB006030B1	6	30
BUSB007012B1	7	12
BUSB007014B1	7	14
BUSB007016B1	7	16
BUSB007018B1	7	18
BUSB007020B1	7	20
BUSB007022B1	7	22
BUSB007024B1	7	24
BUSB007026B1	7	26
BUSB007028B1	7	28
BUSB007030B1	7	30
BUSB008012B1	8	12
BUSB008014B1	8	14
BUSB008016B1	8	16
BUSB008018B1	8	18
BUSB008020B1	8	20
BUSB008022B1	8	22
BUSB008024B1	8	24
BUSB008026B1	8	26
BUSB008028B1	8	28
BUSB008030B1	8	30

P FLEX STENT ONE END CLOSED
WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B2	3	10
BUSB003012B2	3	12
BUSB003014B2	3	14
BUSB003016B2	3	16
BUSB004010B2	4	10
BUSB004012B2	4	12
BUSB004014B2	4	14
BUSB004016B2	4	16
BUSB047012B2	4.7	12
BUSB047014B2	4.7	14
BUSB047016B2	4.7	16
BUSB047018B2	4.7	18
BUSB047020B2	4.7	20
BUSB047022B2	4.7	22
BUSB047024B2	4.7	24
BUSB047026B2	4.7	26
BUSB047028B2	4.7	28
BUSB047030B2	4.7	30
BUSB048012B2	4.8	12
BUSB048014B2	4.8	14
BUSB048016B2	4.8	16
BUSB048018B2	4.8	18
BUSB048020B2	4.8	20
BUSB048022B2	4.8	22
BUSB048024B2	4.8	24
BUSB048026B2	4.8	26
BUSB048028B2	4.8	28
BUSB048030B2	4.8	30
BUSB005012B2	5	12
BUSB005014B2	5	14
BUSB005016B2	5	16
BUSB005018B2	5	18
BUSB005020B2	5	20
BUSB005022B2	5	22
BUSB005024B2	5	24
BUSB005026B2	5	26
BUSB005028B2	5	28
BUSB005030B2	5	30
BUSB006012B2	6	12
BUSB006014B2	6	14
BUSB006016B2	6	16
BUSB006018B2	6	18
BUSB006020B2	6	20
BUSB006022B2	6	22
BUSB006024B2	6	24
BUSB006026B2	6	26
BUSB006028B2	6	28
BUSB006030B2	6	30
BUSB007012B2	7	12
BUSB007014B2	7	14
BUSB007016B2	7	16
BUSB007018B2	7	18
BUSB007020B2	7	20
BUSB007022B2	7	22
BUSB007024B2	7	24
BUSB007026B2	7	26
BUSB007028B2	7	28
BUSB007030B2	7	30
BUSB008012B2	8	12
BUSB008014B2	8	14
BUSB008016B2	8	16
BUSB008018B2	8	18
BUSB008020B2	8	20
BUSB008022B2	8	22
BUSB008024B2	8	24
BUSB008026B2	8	26
BUSB008028B2	8	28
BUSB008030B2	8	30

**P FLEX STENT BOTH ENDS OPEN
WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B3	3	10
BUSB003012B3	3	12
BUSB003014B3	3	14
BUSB003016B3	3	16
BUSB004010B3	4	10
BUSB004012B3	4	12
BUSB004014B3	4	14
BUSB004016B3	4	16
BUSB047012B3	4.7	12
BUSB047014B3	4.7	14
BUSB047016B3	4.7	16
BUSB047018B3	4.7	18
BUSB047020B3	4.7	20
BUSB047022B3	4.7	22
BUSB047024B3	4.7	24
BUSB047026B3	4.7	26
BUSB047028B3	4.7	28
BUSB047030B3	4.7	30
BUSB048012B3	4.8	12
BUSB048014B3	4.8	14
BUSB048016B3	4.8	16
BUSB048018B3	4.8	18
BUSB048020B3	4.8	20
BUSB048022B3	4.8	22
BUSB048024B3	4.8	24
BUSB048026B3	4.8	26
BUSB048028B3	4.8	28
BUSB048030B3	4.8	30
BUSB005012B3	5	12
BUSB005014B3	5	14
BUSB005016B3	5	16
BUSB005018B3	5	18
BUSB005020B3	5	20
BUSB005022B3	5	22
BUSB005024B3	5	24
BUSB005026B3	5	26
BUSB005028B3	5	28
BUSB005030B3	5	30
BUSB006012B3	6	12
BUSB006014B3	6	14
BUSB006016B3	6	16
BUSB006018B3	6	18
BUSB006020B3	6	20
BUSB006022B3	6	22
BUSB006024B3	6	24
BUSB006026B3	6	26
BUSB006028B3	6	28
BUSB006030B3	6	30
BUSB007012B3	7	12
BUSB007014B3	7	14
BUSB007016B3	7	16
BUSB007018B3	7	18
BUSB007020B3	7	20
BUSB007022B3	7	22
BUSB007024B3	7	24
BUSB007026B3	7	26
BUSB007028B3	7	28
BUSB007030B3	7	30
BUSB008012B3	8	12
BUSB008014B3	8	14
BUSB008016B3	8	16
BUSB008018B3	8	18
BUSB008020B3	8	20
BUSB008022B3	8	22
BUSB008024B3	8	24
BUSB008026B3	8	26
BUSB008028B3	8	28
BUSB008030B3	8	30

**P FLEX STENT ONE END CLOSED
WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B4	3	10
BUSB003012B4	3	12
BUSB003014B4	3	14
BUSB003016B4	3	16
BUSB004010B4	4	10
BUSB004012B4	4	12
BUSB004014B4	4	14
BUSB004016B4	4	16
BUSB047012B4	4.7	12
BUSB047014B4	4.7	14
BUSB047016B4	4.7	16
BUSB047018B4	4.7	18
BUSB047020B4	4.7	20
BUSB047022B4	4.7	22
BUSB047024B4	4.7	24
BUSB047026B4	4.7	26
BUSB047028B4	4.7	28
BUSB047030B4	4.7	30
BUSB048012B4	4.8	12
BUSB048014B4	4.8	14
BUSB048016B4	4.8	16
BUSB048018B4	4.8	18
BUSB048020B4	4.8	20
BUSB048022B4	4.8	22
BUSB048024B4	4.8	24
BUSB048026B4	4.8	26
BUSB048028B4	4.8	28
BUSB048030B4	4.8	30
BUSB005012B4	5	12
BUSB005014B4	5	14
BUSB005016B4	5	16
BUSB005018B4	5	18
BUSB005020B4	5	20
BUSB005022B4	5	22
BUSB005024B4	5	24
BUSB005026B4	5	26
BUSB005028B4	5	28
BUSB005030B4	5	30
BUSB006012B4	6	12
BUSB006014B4	6	14
BUSB006016B4	6	16
BUSB006018B4	6	18
BUSB006020B4	6	20
BUSB006022B4	6	22
BUSB006024B4	6	24
BUSB006026B4	6	26
BUSB006028B4	6	28
BUSB006030B4	6	30
BUSB007012B4	7	12
BUSB007014B4	7	14
BUSB007016B4	7	16
BUSB007018B4	7	18
BUSB007020B4	7	20
BUSB007022B4	7	22
BUSB007024B4	7	24
BUSB007026B4	7	26
BUSB007028B4	7	28
BUSB007030B4	7	30
BUSB008012B4	8	12
BUSB008014B4	8	14
BUSB008016B4	8	16
BUSB008018B4	8	18
BUSB008020B4	8	20
BUSB008022B4	8	22
BUSB008024B4	8	24
BUSB008026B4	8	26
BUSB008028B4	8	28
BUSB008030B4	8	30

**P FLEX STENT BOTH ENDS OPEN
BOTH ENDS MULTILOOP WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSB004022B5	4	22-32
BUSB047018B5 BUSB047024B5	4.7 4.7	18 24
BUSB048018B5 BUSB048024B5	4.8 4.8	18 24
BUSB005018B5 BUSB005024B5	5 5	18 24
BUSB006018B5 BUSB006024B5	6 6	18 24
BUSB007018B5 BUSB007024B5	7 7	18 24
BUSB008018B5 BUSB008024B5	8 8	18 24

**P FLEX STENT ONE END CLOSED
BOTH ENDS MULTILOOP WITH STRING**

Product Code	Size (in Fr)	Length (in cm)
BUSB004022B6	4	22-32
BUSB047018B6 BUSB047024B6	4.7 4.7	18 24
BUSB048018B6 BUSB048024B6	4.8 4.8	18 24
BUSB005018B6 BUSB005024B6	5 5	18 24
BUSB006018B6 BUSB006024B6	6 6	18 24
BUSB007018B6 BUSB007024B6	7 7	18 24
BUSB008018B6 BUSB008024B6	8 8	18 24

**P FLEX STENT BOTH ENDS OPEN
BEML WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB047018B7 BUSB047024B7	4.7 4.7	18 24
BUSB048018B7 BUSB048024B7	4.8 4.8	18 24
BUSB005018B7 BUSB005024B7	5 5	18 24
BUSB006018B7 BUSB006024B7	6 6	18 24
BUSB007018B7 BUSB007024B7	7 7	18 24
BUSB008018B7 BUSB008024B7	8 8	18 24

**P FLEX STENT ONE END CLOSED
BEML WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB047018B8 BUSB047024B8	4.7 4.7	18 24
BUSB048018B8 BUSB048024B8	4.8 4.8	18 24
BUSB005018B8 BUSB005024B8	5 5	18 24
BUSB006018B8 BUSB006024B8	6 6	18 24
BUSB007018B8 BUSB007024B8	7 7	18 24
BUSB008018B8 BUSB008024B8	8 8	18 24

P FLEX STENT SET WITH STRING

(URETERAL STENT SET / DOUBLE J STENT SET)

Intended to be implanted for drainage of urine from renal pelvis to urinary bladder.

Key Features


- For indwelling up to 30 days
- Polyurethane based stent has uniform diameter and smooth surface reducing encrustation
- Optimal loop shape and high memory retention reduces the risk of migration
- Optimal radiopacity for better visualization
- Stent with string which helps in repositioning

Set Contains

Stent with string, pusher, guidewire & clamp

Types

- Both Ends Open
- One End Closed
- Both Ends Multiloop
- With or Without Hydrophilic Coating



Standard

Size (in Fr)	Length (in cm)
4.7, 4.8, 5, 6, 7 & 8	12, 14, 16, 18, 20, 22, 24, 26 28 & 30

Multiloop

Size (in Fr)	Length (in cm)
4.7, 4.8, 5, 6, 7 & 8	18-24 & 24-30

Pediatric

Size (in Fr)	Length (in cm)
3 & 4	10, 12, 14 & 16

P FLEX STENT SET BOTH ENDS OPEN
WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B9	3	10
BUSB003012B9	3	12
BUSB003014B9	3	14
BUSB003016B9	3	16
BUSB004010B9	4	10
BUSB004012B9	4	12
BUSB004014B9	4	14
BUSB004016B9	4	16
BUSB047012B9	4.7	12
BUSB047014B9	4.7	14
BUSB047016B9	4.7	16
BUSB047018B9	4.7	18
BUSB047020B9	4.7	20
BUSB047022B9	4.7	22
BUSB047024B9	4.7	24
BUSB047026B9	4.7	26
BUSB047028B9	4.7	28
BUSB047030B9	4.7	30
BUSB048012B9	4.8	12
BUSB048014B9	4.8	14
BUSB048016B9	4.8	16
BUSB048018B9	4.8	18
BUSB048020B9	4.8	20
BUSB048022B9	4.8	22
BUSB048024B9	4.8	24
BUSB048026B9	4.8	26
BUSB048028B9	4.8	28
BUSB048030B9	4.8	30
BUSB005012B9	5	12
BUSB005014B9	5	14
BUSB005016B9	5	16
BUSB005018B9	5	18
BUSB005020B9	5	20
BUSB005022B9	5	22
BUSB005024B9	5	24
BUSB005026B9	5	26
BUSB005028B9	5	28
BUSB005030B9	5	30
BUSB006012B9	6	12
BUSB006014B9	6	14
BUSB006016B9	6	16
BUSB006018B9	6	18
BUSB006020B9	6	20
BUSB006022B9	6	22
BUSB006024B9	6	24
BUSB006026B9	6	26
BUSB006028B9	6	28
BUSB006030B9	6	30
BUSB007012B9	7	12
BUSB007014B9	7	14
BUSB007016B9	7	16
BUSB007018B9	7	18
BUSB007020B9	7	20
BUSB007022B9	7	22
BUSB007024B9	7	24
BUSB007026B9	7	26
BUSB007028B9	7	28
BUSB007030B9	7	30
BUSB008012B9	8	12
BUSB008014B9	8	14
BUSB008016B9	8	16
BUSB008018B9	8	18
BUSB008020B9	8	20
BUSB008022B9	8	22
BUSB008024B9	8	24
BUSB008026B9	8	26
BUSB008028B9	8	28
BUSB008030B9	8	30

P FLEX STENT SET ONE END CLOSED
WITH STRING

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B10	3	10
BUSB003012B10	3	12
BUSB003014B10	3	14
BUSB003016B10	3	16
BUSB004010B10	4	10
BUSB004012B10	4	12
BUSB004014B10	4	14
BUSB004016B10	4	16
BUSB047012B10	4.7	12
BUSB047014B10	4.7	14
BUSB047016B10	4.7	16
BUSB047018B10	4.7	18
BUSB047020B10	4.7	20
BUSB047022B10	4.7	22
BUSB047024B10	4.7	24
BUSB047026B10	4.7	26
BUSB047028B10	4.7	28
BUSB047030B10	4.7	30
BUSB048012B10	4.8	12
BUSB048014B10	4.8	14
BUSB048016B10	4.8	16
BUSB048018B10	4.8	18
BUSB048020B10	4.8	20
BUSB048022B10	4.8	22
BUSB048024B10	4.8	24
BUSB048026B10	4.8	26
BUSB048028B10	4.8	28
BUSB048030B10	4.8	30
BUSB005012B10	5	12
BUSB005014B10	5	14
BUSB005016B10	5	16
BUSB005018B10	5	18
BUSB005020B10	5	20
BUSB005022B10	5	22
BUSB005024B10	5	24
BUSB005026B10	5	26
BUSB005028B10	5	28
BUSB005030B10	5	30
BUSB006012B10	6	12
BUSB006014B10	6	14
BUSB006016B10	6	16
BUSB006018B10	6	18
BUSB006020B10	6	20
BUSB006022B10	6	22
BUSB006024B10	6	24
BUSB006026B10	6	26
BUSB006028B10	6	28
BUSB006030B10	6	30
BUSB007012B10	7	12
BUSB007014B10	7	14
BUSB007016B10	7	16
BUSB007018B10	7	18
BUSB007020B10	7	20
BUSB007022B10	7	22
BUSB007024B10	7	24
BUSB007026B10	7	26
BUSB007028B10	7	28
BUSB007030B10	7	30
BUSB008012B10	8	12
BUSB008014B10	8	14
BUSB008016B10	8	16
BUSB008018B10	8	18
BUSB008020B10	8	20
BUSB008022B10	8	22
BUSB008024B10	8	24
BUSB008026B10	8	26
BUSB008028B10	8	28
BUSB008030B10	8	30

**P FLEX STENT SET BOTH ENDS OPEN
WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B11	3	10
BUSB003012B11	3	12
BUSB003014B11	3	14
BUSB003016B11	3	16
BUSB004010B11	4	10
BUSB004012B11	4	12
BUSB004014B11	4	14
BUSB004016B11	4	16
BUSB047012B11	4.7	12
BUSB047014B11	4.7	14
BUSB047016B11	4.7	16
BUSB047018B11	4.7	18
BUSB047020B11	4.7	20
BUSB047022B11	4.7	22
BUSB047024B11	4.7	24
BUSB047026B11	4.7	26
BUSB047028B11	4.7	28
BUSB047030B11	4.7	30
BUSB048012B11	4.8	12
BUSB048014B11	4.8	14
BUSB048016B11	4.8	16
BUSB048018B11	4.8	18
BUSB048020B11	4.8	20
BUSB048022B11	4.8	22
BUSB048024B11	4.8	24
BUSB048026B11	4.8	26
BUSB048028B11	4.8	28
BUSB048030B11	4.8	30
BUSB005012B11	5	12
BUSB005014B11	5	14
BUSB005016B11	5	16
BUSB005018B11	5	18
BUSB005020B11	5	20
BUSB005022B11	5	22
BUSB005024B11	5	24
BUSB005026B11	5	26
BUSB005028B11	5	28
BUSB005030B11	5	30
BUSB006012B11	6	12
BUSB006014B11	6	14
BUSB006016B11	6	16
BUSB006018B11	6	18
BUSB006020B11	6	20
BUSB006022B11	6	22
BUSB006024B11	6	24
BUSB006026B11	6	26
BUSB006028B11	6	28
BUSB006030B11	6	30
BUSB007012B11	7	12
BUSB007014B11	7	14
BUSB007016B11	7	16
BUSB007018B11	7	18
BUSB007020B11	7	20
BUSB007022B11	7	22
BUSB007024B11	7	24
BUSB007026B11	7	26
BUSB007028B11	7	28
BUSB007030B11	7	30
BUSB008012B11	8	12
BUSB008014B11	8	14
BUSB008016B11	8	16
BUSB008018B11	8	18
BUSB008020B11	8	20
BUSB008022B11	8	22
BUSB008024B11	8	24
BUSB008026B11	8	26
BUSB008028B11	8	28
BUSB008030B11	8	30

**P FLEX STENT SET ONE END CLOSED
WITH STRING & HYDROPHILIC COATING**

Product Code	Size (in Fr)	Length (in cm)
BUSB003010B12	3	10
BUSB003012B12	3	12
BUSB003014B12	3	14
BUSB003016B12	3	16
BUSB004010B12	4	10
BUSB004012B12	4	12
BUSB004014B12	4	14
BUSB004016B12	4	16
BUSB047012B12	4.7	12
BUSB047014B12	4.7	14
BUSB047016B12	4.7	16
BUSB047018B12	4.7	18
BUSB047020B12	4.7	20
BUSB047022B12	4.7	22
BUSB047024B12	4.7	24
BUSB047026B12	4.7	26
BUSB047028B12	4.7	28
BUSB047030B12	4.7	30
BUSB048012B12	4.8	12
BUSB048014B12	4.8	14
BUSB048016B12	4.8	16
BUSB048018B12	4.8	18
BUSB048020B12	4.8	20
BUSB048022B12	4.8	22
BUSB048024B12	4.8	24
BUSB048026B12	4.8	26
BUSB048028B12	4.8	28
BUSB048030B12	4.8	30
BUSB005012B12	5	12
BUSB005014B12	5	14
BUSB005016B12	5	16
BUSB005018B12	5	18
BUSB005020B12	5	20
BUSB005022B12	5	22
BUSB005024B12	5	24
BUSB005026B12	5	26
BUSB005028B12	5	28
BUSB005030B12	5	30
BUSB006012B12	6	12
BUSB006014B12	6	14
BUSB006016B12	6	16
BUSB006018B12	6	18
BUSB006020B12	6	20
BUSB006022B12	6	22
BUSB006024B12	6	24
BUSB006026B12	6	26
BUSB006028B12	6	28
BUSB006030B12	6	30
BUSB007012B12	7	12
BUSB007014B12	7	14
BUSB007016B12	7	16
BUSB007018B12	7	18
BUSB007020B12	7	20
BUSB007022B12	7	22
BUSB007024B12	7	24
BUSB007026B12	7	26
BUSB007028B12	7	28
BUSB007030B12	7	30
BUSB008012B12	8	12
BUSB008014B12	8	14
BUSB008016B12	8	16
BUSB008018B12	8	18
BUSB008020B12	8	20
BUSB008022B12	8	22
BUSB008024B12	8	24
BUSB008026B12	8	26
BUSB008028B12	8	28
BUSB008030B12	8	30

P FLEX STENT SET BOTH ENDS OPEN BEML WITH STRING			P FLEX STENT SET ONE END CLOSED BEML WITH STRING		
Product Code	Size (in Fr)	Length (in cm)	Product Code	Size (in Fr)	Length (in cm)
BUSB004022B13	4	22-32	BUSB004022B14	4	22-32
BUSB047018B13	4.7	18	BUSB047018B14	4.7	18
BUSB047024B13	4.7	24	BUSB047024B14	4.7	24
BUSB048018B13	4.8	18	BUSB048018B14	4.8	18
BUSB048024B13	4.8	24	BUSB048024B14	4.8	24
BUSB005018B13	5	18	BUSB005018B14	5	18
BUSB005024B13	5	24	BUSB005024B14	5	24
BUSB006018B13	6	18	BUSB006018B14	6	18
BUSB006024B13	6	24	BUSB006024B14	6	24
BUSB007018B13	7	18	BUSB007018B14	7	18
BUSB007024B13	7	24	BUSB007024B14	7	24
BUSB008018B13	8	18	BUSB008018B14	8	18
BUSB008024B13	8	24	BUSB008024B14	8	24

P FLEX STENT SET BOTH ENDS OPEN BEML WITH STRING & HYDROPHILIC COATING			P FLEX STENT SET ONE END CLOSED BEML WITH STRING & HYDROPHILIC COATING		
Product Code	Size (in Fr)	Length (in cm)	Product Code	Size (in Fr)	Length (in cm)
BUSB047018B15	4.7	18	BUSB047018B16	4.7	18
BUSB047024B15	4.7	24	BUSB047024B16	4.7	24
BUSB048018B15	4.8	18	BUSB048018B16	4.8	18
BUSB048024B15	4.8	24	BUSB048024B16	4.8	24
BUSB005018B15	5	18	BUSB005018B16	5	18
BUSB005024B15	5	24	BUSB005024B16	5	24
BUSB006018B15	6	18	BUSB006018B16	6	18
BUSB006024B15	6	24	BUSB006024B16	6	24
BUSB007018B15	7	18	BUSB007018B16	7	18
BUSB007024B15	7	24	BUSB007024B16	7	24
BUSB008018B15	8	18	BUSB008018B16	8	18
BUSB008024B15	8	24	BUSB008024B16	8	24

ENDOPYELOTOMY STENT

Used for internal drainage from uretero pelvic junction to the bladder following incision of a stricture.

Key Features

- Dual diameter design
- Smooth transition from larger diameter to smaller diameter
- Maintain ureteral patency post Endopyelotomy

Also available in Set, which Includes


Stent, pusher, guidewire & clamp

Size (in Fr)

6/12 & 7/14

Length (in cm)

24 & 26



ENDOPYELOTOMY STENT

Product Code	Size (in Fr)	Length (in cm)
BEPB6I2O24B1	6/12	24
BEPB6I2O26B1	6/12	26
BEPB7I4O24B1	7/14	24
BEPB7I4O26B1	7/14	26

ENDOPYELOTOMY STENT SET

Product Code	Size (in Fr)	Length (in cm)
BEPB6I2O24B2	6/12	24
BEPB6I2O26B2	6/12	26
BEPB7I4O24B2	7/14	24
BEPB7I4O26B2	7/14	26

URETHRAL STENT

Used for stenting the urethra after Hypospedias or Epispedias repair and to provide post-operative drainage of the bladder.

Key Features


- Smooth surface for ease of introduction
- Radiopaque for better visualization
- Hydrophilic coating for easy insertion

Size (in Fr)

5, 6, 7, 8, 10 & 12

Length (in cm)

25



URETHRAL STENT OPEN END

Product Code	Size (in Fr)	Length (in cm)
BHSB005025B1	5	25
BHSB006025B1	6	25
BHSB007025B1	7	25
BHSB008025B1	8	25
BHSB010025B1	10	25
BHSB012025B1	12	25

URETHRAL STENT CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BHSB005025B2	5	25
BHSB006025B2	6	25
BHSB007025B2	7	25
BHSB008025B2	8	25
BHSB010025B2	10	25
BHSB012025B2	12	25

URETERAL ACCESS SHEATH

- Praveza[®] Ureteral Access Sheath

UROLOGY

PRAVEZA[®] URETERAL ACCESS SHEATH

Used to establish a conduit during urological procedures facilitating the passage for Ureteroscope / Endoscope and other surgical instruments into the urinary tract (RIRS).

Key Features

- Kink resistant – reinforced with a stainless steel coil designed to provide stability and strength through tortuous anatomy
- Tapered tip dilator with hydrophilic coating provides easy access
- Transparent tip with platinum ring helps in better visualization and placement of sheath under fluoroscopy
- Sizes available for all types of flexible scopes (Storz, Wolf & Single use flexible scopes)
- Seamless transition between dilator and sheath for easy and atraumatic insertion

Size (in Fr)

ID	9.5	10	11.5	12
OD	11.5	12	13.5	14

Length (in cm)

28, 35, 45 & 55



PRAVEZA[®] URETERAL ACCESS SHEATH

Product Code	Size (in Fr)	Length (in cm)
BUAB095028B1	9.5	28
BUAB095035B1	9.5	35
BUAB095045B1	9.5	45
BUAB095055B1	9.5	55
BUAB010028B1	10	28
BUAB010035B1	10	35
BUAB010045B1	10	45
BUAB010055B1	10	55
BUAB115028B1	11.5	28
BUAB115035B1	11.5	35
BUAB115045B1	11.5	45
BUAB115055B1	11.5	55
BUAB012028B1	12	28
BUAB012035B1	12	35
BUAB012045B1	12	45
BUAB012055B1	12	55

CATHETERS

- Ureteral Catheter
- Mono J Catheter
- PCN Catheter
- PCN Catheter with Locking Mechanism
- PCN Catheter with Needle
(Nephrostomy Drainage Set Direct Puncture)
- Malecot Catheter
- Re-entry Malecot Drainage Catheter
- Suprapubic Catheter (SPC)
- CI Catheter

UROLOGY

URETERAL CATHETER

- Used for catheterization of the urinary tract during intraoperative procedures
- Used for injection of contrast media

Key Features

- Optimal radiopacity for better visualization
- With 1 cm marking which helps in catheter position, measurement of ureter length
- Detachable hub for scope insertion

Types


- Open End
- Closed End
- Cone Tip
- Angled Tip

Size (in Fr)

3, 4, 5, 6, & 7

Length (in cm)

70



URETERAL CATHETER STRAIGHT OPEN END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B1	3	70
BUCP004070B1	4	70
BUCP005070B1	5	70
BUCP006070B1	6	70
BUCP007070B1	7	70

URETERAL CATHETER STRAIGHT CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B2	3	70
BUCP004070B2	4	70
BUCP005070B2	5	70
BUCP006070B2	6	70
BUCP007070B2	7	70

URETERAL CATHETER CONE TIP OPEN END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B3	3	70
BUCP004070B3	4	70
BUCP005070B3	5	70
BUCP006070B3	6	70
BUCP007070B3	7	70

URETERAL CATHETER CONE TIP CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B4	3	70
BUCP004070B4	4	70
BUCP005070B4	5	70
BUCP006070B4	6	70
BUCP007070B4	7	70

URETERAL CATHETER ANGLED TIP OPEN END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B5	3	70
BUCP004070B5	4	70
BUCP005070B5	5	70
BUCP006070B5	6	70
BUCP007070B5	7	70

URETERAL CATHETER ANGLED TIP CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BUCP003070B6	3	70
BUCP004070B6	4	70
BUCP005070B6	5	70
BUCP006070B6	6	70
BUCP007070B6	7	70

MONO J CATHETER

Used for temporary drainage and irrigation of urine from the kidney (Pigtail ureteral catheter).

Key Features

- Pigtail loop helps in retention of distal tip in the renal pelvis
- Holes in shaft helps in better drainage
- Used to inject contrast
- Used for short term external drainage to avoid stenting

Also available in Set, which Includes

Mono J catheter, guidewire & urine bag connector

Type


- Open End
- Closed End

Size (in Fr)

4.7, 5, 6, 7 & 8

Length (in cm)

70



MONO J CATHETER OPEN END

Product Code	Size (in Fr)	Length (in cm)
BMJBO47070B1	4.7	70
BMJBO05070B1	5	70
BMJBO06070B1	6	70
BMJBO07070B1	7	70
BMJBO08070B1	8	70

MONO J CATHETER CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BMJBO47070B2	4.7	70
BMJBO05070B2	5	70
BMJBO06070B2	6	70
BMJBO07070B2	7	70
BMJBO08070B2	8	70

MONO J CATHETER SET OPEN END

Product Code	Size (in Fr)	Length (in cm)
BMJBO47070B3	4.7	70
BMJBO05070B3	5	70
BMJBO06070B3	6	70
BMJBO07070B3	7	70
BMJBO08070B3	8	70

MONO J CATHETER SET CLOSED END

Product Code	Size (in Fr)	Length (in cm)
BMJBO47070B4	4.7	70
BMJBO05070B4	5	70
BMJBO06070B4	6	70
BMJBO07070B4	7	70
BMJBO08070B4	8	70

PCN CATHETER

PCN catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion

Package Includes

PCN catheter with metal & plastic straightener and urine bag connector

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

30



PCN CATHETER			PCN CATHETER (BALLOON TYPE)	
Product Code	Size (in Fr)	Length (in cm)	Product Code	Size (in Fr)
BPCP006030B1	6	30	BPCP010000B8	10
BPCP007030B1	7	30	BPCP012000B8	12
BPCP008030B1	8	30	BPCP014000B8	14
BPCP010030B1	10	30	BPCP016000B8	16
BPCP012030B1	12	30	BPCP018000B8	18
BPCP014030B1	14	30	BPCP020000B8	20
BPCP016030B1	16	30		

PCN CATHETER WITH LOCKING MECHANISM

PCN catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine. The locking mechanism prevents accidental slippage of the catheter.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in renal pelvis
- The locking mechanism prevents accidental slippage of the catheter
- Eye holes for better drainage
- Hydrophilic coating for easy insertion

Package Includes

PCN catheter with locking mechanism, metal & plastic straightener and urine bag connector

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

35



Note: Once the catheter is deployed, pull the string, wind and lock to the hub of the catheter.

PCN CATHETER WITH LOCKING MECHANISM

Product Code	Size (in Fr)	Length (in cm)
BPCP006035B2	6	35
BPCP007035B2	7	35
BPCP008035B2	8	35
BPCP010035B2	10	35
BPCP012035B2	12	35
BPCP014035B2	14	35
BPCP016035B2	16	35

PCN CATHETER WITH NEEDLE

(NEPHROSTOMY DRAINAGE SET DIRECT PUNCTURE)

Used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage by direct puncture.

Key Features

- The pigtail tip configuration helps in the retention of the distal tip in the renal pelvis
- Eye holes for better drainage
- Hydrophilic coating for easy insertion
- Trocar design and catheter design prevents buckling of catheter during the procedure
- Available with or without locking mechanism

Package Includes


PCN catheter with metal straightner, metal trocar, urine bag connector and scalpel blade

PCN CATHETER WITH NEEDLE

Size (in Fr)	Length (in cm)
7, 8, 9, 10, 12, 14 16 & 18	30

PCN CATHETER WITH NEEDLE & LOCKING MECHANISM

Size (in Fr)	Length (in cm)
7, 8, 10, 12, 14 & 16	35



PCN CATHETER WITH NEEDLE

Product Code	Size (in Fr)	Length (in cm)
BPCP007030B5	7	30
BPCP008030B5	8	30
BPCP009030B5	9	30
BPCP010030B5	10	30
BPCP012030B5	12	30
BPCP014030B5	14	30
BPCP016030B5	16	30
BPCP018030B5	18	30

PCN CATHETER WITH NEEDLE & LOCKING MECHANISM

Product Code	Size (in Fr)	Length (in cm)
BPCP007035B6	7	35
BPCP008035B6	8	35
BPCP010035B6	10	35
BPCP012035B6	12	35
BPCP014035B6	14	35
BPCP016035B6	16	35

MALECOT CATHETER

Malecot catheter is intended to be inserted percutaneously into the renal pelvis to perform external drainage of urine.

Key Features

- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion

Package Includes

Malecot catheter with plastic & metal straightner and urine bag connector

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

30



MALECOT CATHETER

Product Code	Size (in Fr)	Length (in cm)
BMCP008030B1	8	30
BMCP010030B1	10	30
BMCP012030B1	12	30
BMCP014030B1	14	30
BMCP016030B1	16	30
BMCP018030B1	18	30
BMCP020030B1	20	30
BMCP022030B1	22	30

RE-ENTRY MALECOT DRAINAGE CATHETER

Used to be inserted percutaneously into the renal pelvis upto the ureter to perform external drainage of urine from the renal pelvis and internal drainage into the bladder.

Key Features

- The Malecot flower design is ideally suitable for the drainage of thick viscous fluids
- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Used as a nephrostomy catheter and also as ureteral stent

Package Includes

Re-Entry Malecot Drainage Catheter with plastic straightener & urine bag connector

Size (in Fr)

10/6, 12/6, 14/6, 16/8 & 18/8

Length (in cm)

45



RE-ENTRY MALECOT DRAINAGE CATHETER

Product Code	Size (in Fr)	Length (in cm)
BRMN610045B1	10/6	45
BRMN612045B1	12/6	45
BRMN614045B1	14/6	45
BRMN816045B1	16/8	45
BRMN818045B1	18/8	45

SUPRAPUBIC CATHETER (SPC)

Suprapubic catheter is intended to be inserted percutaneously into the urinary bladder and perform external drainage of urine from the urinary bladder.

Key Features

- Integral flower design matches with the catheter size when closed and cause less trauma during insertion
- Trocar and catheter designed to prevent buckling of catheter during the procedure

Package Includes


Suprapubic catheter with metal straightner, metal trocar and urine bag connector

Size (in Fr)

8, 10, 12, 14, 16 & 18

Length (in cm)

25



SUPRAPUBIC CATHETER

Product Code	Size (in Fr)	Length (in cm)
BSCP008025B1	8	25
BSCP010025B1	10	25
BSCP012025B1	12	25
BSCP014025B1	14	25
BSCP016025B1	16	25
BSCP018025B1	18	25

CI CATHETER

Used to drain urine from bladder.

Key Features

- Easy advancement minimizing the potential trauma
- Firm enough for easy insertion, yet gentle on soft urethral tissue
- Eye holes for better drainage

CI Catheter Female

Size (in Fr)	Length (in cm)
10, 12, 14, 16 & 18	15

CI Catheter Male

Size (in Fr)	Length (in cm)
10, 12, 14, 16 & 18	35



CI CATHETER FEMALE

Product Code	Size (in Fr)	Length (in cm)
BCI0010015B0	10	15
BCI0012015B0	12	15
BCI0014015B0	14	15
BCI0016015B0	16	15
BCI0018015B0	18	15

CI CATHETER MALE

Product Code	Size (in Fr)	Length (in cm)
BCI0010035B0	10	35
BCI0012035B0	12	35
BCI0014035B0	14	35
BCI0016035B0	16	35
BCI0018035B0	18	35

DRAINAGE SETS

- PCN Set - Percutaneous Nephrostomy Drainage Set "Pigtail"
- PCN Set - Pigtail with Locking Mechanism
- PCN Set - Percutaneous Nephrostomy Drainage Set "Malecot"
- Suprapubic Drainage Catheter Set

UROLOGY

PCN SET (PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL")

PCN Set is used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Contains

IP Needle | Guidewire J Tip PTFE Coated
| Fascial Dilators | Urine Bag Connector
PCN Catheter with Hydrophilic Coating
Scalpel Blade | Catheter Fixing Disc

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

30



PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL"

Product Code	Size (in Fr)	Length (in cm)
BPCP006030B3	6	30
BPCP007030B3	7	30
BPCP008030B3	8	30
BPCP010030B3	10	30
BPCP012030B3	12	30
BPCP014030B3	14	30
BPCP016030B3	16	30

PCN SET - PIGTAIL WITH LOCKING MECHANISM

(PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL")

PCN Set is used for percutaneous placement of a pigtail catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Contains

IP Needle | Guidewire J Tip PTFE Coated
| Fascial Dilators | Urine Bag Connector
PCN Catheter with Locking Mechanism
& Hydrophilic Coating | Scalpel Blade
Catheter Fixing Disc

Size (in Fr)

6, 7, 8, 10, 12, 14 & 16

Length (in cm)

35



PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "PIGTAIL" WITH LOCKING MECHANISM

Product Code	Size (in Fr)	Length (in cm)
BPCP006035B4	6	35
BPCP007035B4	7	35
BPCP008035B4	8	35
BPCP010035B4	10	35
BPCP012035B4	12	35
BPCP014035B4	14	35
BPCP016035B4	16	35

PCN SET (PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "MALECOT")

PCN Set is used for percutaneous placement of a Malecot catheter in the renal pelvis for nephrostomy drainage.

Advantages

Kit provides entire devices for performing Nephrostomy procedure

Set Contains

IP Needle | Guidewire J Tip PTFE Coated
| Fascial Dilators | Urine Bag Connector
Malecot Catheter with Hydrophilic
Coating | Scalpel Blade | Catheter Fixing
Disc

Size (in Fr)

8, 10, 12, 14, 16, 18, 20 & 22

Length (in cm)

30



PERCUTANEOUS NEPHROSTOMY DRAINAGE SET "MALECOT"

Product Code	Size (in Fr)	Length (in cm)
BMCP008030B2	8	30
BMCP010030B2	10	30
BMCP012030B2	12	30
BMCP014030B2	14	30
BMCP016030B2	16	30
BMCP018030B2	18	30
BMCP020030B2	20	30
BMCP022030B2	22	30

SUPRAPUBIC DRAINAGE CATHETER SET (SDS)

Suprapubic drainage catheter set is intended to create suprapubic access to the bladder and place a catheter for bladder drainage.

Advantages

- Silicone foley catheter for patient comfort and safety
- Specially designed trocar with needle tip enhances smooth puncture
- Peel away sheath forms a conduit for catheter placement and can be peeled easily to remove from the catheter
- Seamless transition between trocar and peel away sheath for easy and atraumatic puncture

Set Contains

Foley catheter, trocar with peel-away sheath and scalpel blade

Size (in Fr)

10, 12, 14, 16 & 18



SUPRAPUBIC DRAINAGE CATHETER SET - SILICONE

Product Code	Size (in Fr)
BSD0010025B1	10
BSD0012025B1	12
BSD0014025B1	14
BSD0016025B1	16
BSD0018025B1	18

STONE BASKETS

- N Quadra™ - Nitinol Stone Basket With Tip
- N Quadra™ Z - Nitinol Stone Basket Zero Tip
- N Quadra™ Helical - Nitinol Stone Basket Helical With Tip
- N Quadra™ Helical Z - Nitinol Stone Basket Helical Zero Tip
- N Quadra™ Flexible - Nitinol Stone Basket Flexible
- Delta Basket
- PCNL Basket

UROLOGY

N QUADRA™ - NITINOL STONE BASKET WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly
- Atraumatic tip

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ - NITINOL STONE BASKET WITH TIP

Product Code	Size (in Fr)	Length (in cm)
BSBN003070BT1	3	70
BSBN003090BT1	3	90
BSBN003120BT1	3	120
BSBN004070BT1	4	70
BSBN004090BT1	4	90
BSBN004120BT1	4	120

N QUADRA™ Z - NITINOL STONE BASKET ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ Z - NITINOL STONE BASKET ZERO TIP

Product Code	Size (in Fr)	Length (in cm)
BSBN003070BTO	3	70
BSBN003090BTO	3	90
BSBN003120BTO	3	120
BSBN004070BTO	4	70
BSBN004090BTO	4	90
BSBN004120BTO	4	120

N QUADRA™ HELICAL - NITINOL STONE BASKET HELICAL WITH TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ HELICAL - NITINOL STONE BASKET HELICAL WITH TIP

Product Code	Size (in Fr)	Length (in cm)
BSBN003070BHI	3	70
BSBN003090BHI	3	90
BSBN003120BHI	3	120
BSBN004070BHI	4	70
BSBN004090BHI	4	90
BSBN004120BHI	4	120

N QUADRA™ HELICAL Z - NITINOL STONE BASKET HELICAL ZERO TIP

Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from the Ureter or the Kidney.

Key Features

- Zero tip design reduces chances of perforation and mucosal damage, which also avoids stone migration during the procedure
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- Helical design helps in capturing small stones and fragments effectively
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

3 & 4

Length (in cm)

70, 90 & 120



N QUADRA™ HELICAL Z - NITINOL STONE BASKET HELICAL ZERO TIP

Product Code	Size (in Fr)	Length (in cm)
BSBN003070BHO	3	70
BSBN003090BHO	3	90
BSBN003120BHO	3	120
BSBN004070BHO	4	70
BSBN004090BHO	4	90
BSBN004120BHO	4	120

N QUADRA™ FLEXIBLE - NITINOL STONE BASKET FLEXIBLE

Flexible Stone Basket is intended to be used in an endoscopic procedure for removal of urinary stones or foreign objects from calyx.

Key Features

- Flexible design shaft which can angulate up to 270 degrees as per the flexible scope direction during RIRS Procedure
- Kink resistant braided shaft provides durability, making removal of several fragments possible
- Tipless design helps in easy capturing of stones by allowing the basket to open at the base of calyx and encircle the stone
- Nitinol wire construction with high strength and high memory helps the basket to retain its shape during extreme torsion
- In case of stone impregnation in the basket during the procedure, ergonomic handle permits full disassembly

Size (in Fr)

1.8, 1.9, 2, 2.2 & 2.4

Length (in cm)

120



N QUADRA™ FLEXIBLE - NITINOL STONE BASKET FLEXIBLE

Product Code	Size (in Fr)	Length (in cm)
BFBB018120BTO	1.8	120
BFBB019120BTO	1.9	120
BFBB002120BTO	2	120
BFBB022120BTO	2.2	120
BFBB024120BTO	2.4	120

DELTA BASKET

Steerable Basket for rigid nephroscope for stone extraction from calices which are inaccessible due to acute angle.

Design Features

- Can be angled inside kidney to reach inner calyx
- Can be used for straight path stone manipulation and also can be angled both directions
- Unique basket shape which prevents stone migration

Procedural Features

- Prevents injury to the infundibulum and reduce bleeding
- Avoids the need for multiple tracts
- Can be used without the assistance of flexible nephroscope
- Helps to remove big blood clots

Size (in Fr)

7

Length (in cm)

40



DELTA BASKET

Product Code	Size (in Fr)	Length (in cm)
BDBR007040B0	7	40

PCNL BASKET

PCNL Basket is intended to be used with rigid nephroscope for stone extraction from calices.

Design Features

- Easy to handle - single handed operation
- Can be used for straight path stone manipulation
- Unique basket shape which prevents stone migration

Size (in Fr)

7 & 9

Length (in cm)

40



PCNL BASKET

Product Code	Size (in Fr)	Length (in cm)
BPBR007040B0	7	40
BPBR009040B0	9	40

DILATORS

- Fascial Dilator Set
- Ureteral Dilator Set
- Single Step™ Ureteral Dilator
- Urethral Dilator Set
- Urethral Dilator Set - S Curved
- Urethral Self Dilator
- Amplatz Sheath
- Amplatz Renal Dilator
- Amplatz Renal Dilation Set - S
- Amplatz Renal Dilation Set - F
- Mini Perc Renal Dilation Set
- Ureteral Balloon Dilator
- Nephrostomy Balloon Dilator With Amplatz Sheath
- Screw Dilator

UROLOGY

FASCIAL DILATOR SET

Used for fascial dilation, intended for stretching or enlarging the fascial tissue covering a cavity, tract or opening prior to a diagnostic or interventional procedure.

Key Features


- Gradual tapered end & smooth surface for atraumatic entry
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

6, 8, 10, 12, 14, 16, 18 & 20

Length (in cm)

20 & 22



FASCIAL DILATOR - SINGLE

Product Code	Size (in Fr)	Length (in cm)
BFDN006020BS	6	20
BFDN008020BS	8	20
BFDN010020BS	10	20
BFDN012020BS	12	20
BFDN014020BS	14	20
BFDN016020BS	16	20
BFDN018020BS	18	20
BFDN020020BS	20	20
BFDN006022BS	6	22
BFDN008022BS	8	22
BFDN010022BS	10	22
BFDN012022BS	12	22
BFDN014022BS	14	22
BFDN016022BS	16	22
BFDN018022BS	18	22
BFDN020022BS	20	22

FASCIAL DILATOR SET

Product Code	Size (in Fr)	Length (in cm)
BFDN616020BSE	6-16	20
BFDN620020BSE	6-20	20
BFDN616022BSE	6-16	22
BFDN620022BSE	6-20	22

URETERAL DILATOR SET

Used for dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

Key Features


- Tapered end & smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

6, 8, 10, 12, 14 & 16

Length (in cm)

70



URETERAL DILATOR - SINGLE

Product Code	Size (in Fr)	Length (in cm)
BUDN006070BS	6	70
BUDN008070BS	8	70
BUDN010070BS	10	70
BUDN012070BS	12	70
BUDN014070BS	14	70
BUDN016070BS	16	70

URETERAL DILATOR SET

Product Code	Size (in Fr)	Length (in cm)
BUDN616070BSE	6-16	70

SINGLE STEP™ URETERAL DILATOR

Used for one step dilation of the ureter through natural orifice prior to a diagnostic or interventional procedure.

Key Features

- Long tapered tip avoids use of sequential dilation
- Smooth surface for atraumatic entry during dilation
- Also available with hydrophilic coating for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

10/6, 12/6, 14/6, 16/6

Length (in cm)

70



SINGLE STEP™ URETERAL DILATOR		
Product Code	Size (in Fr)	Length (in cm)
BNDP610070BO	10/6	70
BNDP612070BO	12/6	70
BNDP614070BO	14/6	70
BNDP616070BO	16/6	70

URETHRAL DILATOR SET

Used for dilation of the Urethra.

Key Features

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization
- To be used for sequential dilation
- Side hole for drainage

Type


Straight

Size (in Fr)

8-24 (8, 10, 12, 14, 16,18, 20, 22 & 24)
8-30 (8, 10, 12, 14, 16,18, 20, 22, 24, 26, 28 & 30)

Length (in cm)

40



URETHRAL DILATOR SET		
Product Code	Size (in Fr)	Length (in cm)
BUD0824040BSE BUD0830040BSE	8-24 8-30	40 40

URETHRAL DILATOR SET - S CURVED

Used for dilation of the Male Urethral strictures.

Key Features

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization
- To be used for sequential dilation
- Side hole for drainage
- S curve supports natural curvature of male urethra

Type


S Curved

Size (in Fr)

8-24 (8, 10, 12, 14, 16,18, 20, 22 & 24)
8-30 (8, 10, 12, 14, 16,18, 20, 22, 24, 26, 28 & 30)

Length (in cm)

40



URETHRAL DILATOR SET - S CURVED		
Product Code	Size (in Fr)	Length (in cm)
BUD0824040BSS	8-24	40
BUD0830040BSS	8-30	40

URETHRAL SELF DILATOR

Used for self-dilation of Urethral strictures.

Key Features

- Tapered end for atraumatic entry during dilation
- Smooth surface for ease of introduction

Types

- Straight With Hole
- Curved With Hole

Size (in Fr)

6, 8, 10, 12, 14 & 16

Length (in cm)

20 & 30



URETHRAL SELF DILATOR - STRAIGHT WITH HOLE

Product Code	Size (in Fr)	Length (in cm)
BSD0610030B1	10	30
BSD0612030B1	12	30
BSD0614030B1	14	30
BSD0616030B1	16	30

URETHRAL SELF DILATOR - CURVED WITH HOLE

Product Code	Size (in Fr)	Length (in cm)
BSD0610030B3	10	30
BSD0612030B3	12	30
BSD0614030B3	14	30
BSD0616030B3	16	30

AMPLATZ SHEATH

Intended for percutaneous dilation of tissue surrounding a dilator previously placed. It is used as a working channel to maintain a previously established nephrostomy tract for introduction of instruments and catheters during urological procedures.

Key Features

- Optimal bevel tip for atraumatic entry
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization

Size (in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

17 & 20



AMPLATZ SHEATH

Product Code	Size (in Fr)	Length (in cm)
BAST012017B1	12	17
BAST014017B1	14	17
BAST016017B1	16	17
BAST018017B1	18	17
BAST020017B1	20	17
BAST022017B1	22	17
BAST024017B1	24	17
BAST026017B1	26	17
BAST028017B1	28	17
BAST030017B1	30	17
BAST012020B1	12	20
BAST014020B1	14	20
BAST016020B1	16	20
BAST018020B1	18	20
BAST020020B1	20	20
BAST022020B1	22	20
BAST024020B1	24	20
BAST026020B1	26	20
BAST028020B1	28	20
BAST030020B1	30	20

AMPLATZ RENAL DILATOR

(RENAL DILATOR WITH AMPLATZ SHEATH)

Used for dilation of a percutaneous tract and setup of a working channel into the kidney for PCNL procedures by placing a suitable amplatz sheath.

Key Features

- Gradual taper designed for atraumatic placement and maneuverability
- Smooth surface for ease of introduction
- Optimal radiopacity for better visualization

Dilator Size (OD of Dilator in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

30

Amplatz Sheath Size (ID of Sheath in Fr)

12, 14, 16, 18, 20, 22, 24, 26, 28 & 30

Length (in cm)

17



AMPLATZ RENAL DILATOR

Product Code	Size (in Fr)	Length (in cm)
BRD0012030BS	12	30
BRD0014030BS	14	30
BRD0016030BS	16	30
BRD0018030BS	18	30
BRD0020030BS	20	30
BRD0022030BS	22	30
BRD0024030BS	24	30
BRD0026030BS	26	30
BRD0028030BS	28	30
BRD0030030BS	30	30

AMPLATZ RENAL DILATION SET - S

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.


Set Includes

- 8 Fr guiding catheter
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr – 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr – 30 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 24, 26, 28 & 30 Fr of length 17 cm
- Packed in tray

AMPLATZ RENAL DILATION SET - S

Product Code

BAR0630000BSE



AMPLATZ RENAL DILATION SET - F

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.


Set Includes

- Scalpel blade
- 2 part 18 G x 15 cm IP needle with echogenic tip
- 8 Fr guiding catheter
- 0.035" PTFE extra stiff J tip guidewire
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr – 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr – 30 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 24, 26, 28 & 30 Fr of length 17 cm
- Packed in tray

AMPLATZ RENAL DILATION SET - F

Product Code

BAR0630000B53



MINI PERC RENAL DILATION SET

Intended for stretching or enlarging the fascial tissue prior to a nephrostomy procedure and to introduce the amplatz sheath which is intended for percutaneous dilation of tissue surrounding a dilator previously placed and used as a working channel during urological procedures.


Set Includes

- Scalpel blade
- Trocar Needle with sheath 6 Fr x 15 cm
- 8 Fr guiding catheter with compatible 0.035" PTFE extra stiff J tip guidewire
- 10 Fr radiopaque dilation sheath
- Radiopaque fascial dilators of sizes from 6 Fr - 10 Fr of length 20 cm
- Radiopaque renal dilators of sizes from 12 Fr - 22 Fr of length 30 cm
- Teflon amplatz sheath for renal dilators of sizes 14, 16, 18, 20 & 22 Fr of length 17 cm
- Packed in tray

MINI PERC RENAL DILATION SET

Product Code

BAM0622000BSE



URETERAL BALLOON DILATOR


Used for dilation of narrow ureteral segments.

Product Features

- Controlled inflation pressure allows dilating force at a stricture or ureteral narrowing
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy

Product Specifications

- Shaft Size: 5 Fr
- Catheter Length: 75 cm
- Balloon Length: 4 cm
- Balloon Expansion Diameter: 5 mm
- Nominal Pressure: 7 atm
- Rate Burst Pressure: 14 atm
- Use Guidewire: 0.035 inch



URETERAL BALLOON DILATOR			URETERAL BALLOON DILATOR WITH INFLATION DEVICE		
Product Code	Size (in Fr)	Length (in cm)	Product Code	Size (in Fr)	Length (in cm)
BUB0005075B0	5	75	BUB0005075B1	5	75

NEPHROSTOMY BALLOON DILATOR WITH AMPLATZ SHEATH


Used for dilation of percutaneous tract for nephrostomy procedure.

Product Features

- Controlled inflation pressure allows dilating force at a nephrostomy tract
- Radiopaque marker bands enhance visualization of balloon position under fluoroscopy

Product Specifications

- Shaft Size: 5 Fr
- Catheter Length: 55 cm
- Balloon Length: 12 cm
- Balloon Expansion Diameter: 10 mm
- Nominal Pressure: 2 atm
- Rate Burst Pressure: 11 atm
- Use Guidewire: 0.035 inch



NEPHROSTOMY BALLOON DILATOR

Product Code	Size (in Fr)	Length (in cm)
BNB0005055B0	5	55

NEPHROSTOMY BALLOON DILATOR WITH AMPLATZ SHEATH

Product Code	Size (in Fr)	Length (in cm)
BNB0005055B1	5	55

NEPHROSTOMY BALLOON DILATOR WITH AMPLATZ SHEATH & INFLATION DEVICE

Product Code	Size (in Fr)	Length (in cm)
BNB0005055B2	5	55

SCREW DILATOR

Used for fascial dilation in one step.

Key Features

- Radiopaque for better visualization
- Screw tip for one step dilation of fascia during mini PCNL procedure
- Also available with compatible amplatz sheath
- Avoids sequential dilatation, reducing procedural time

Size (in Fr)

12, 14 & 16

Length (in cm)

25



SCREW DILATOR			SCREW DILATOR WITH AMPLATZ SHEATH		
Product Code	Size (in Fr)	Length (in cm)	Product Code	Size (in Fr)	Length (in cm)
BOS0012025B0	12	25	BOS0012025BAS	12	25
BOS0014025B0	14	25	BOS0014025BAS	14	25
BOS0016025B0	16	25	BOS0016025BAS	16	25

NEEDLES

- Initial Puncture Needle - Two Part
- Sure Puncture™ Needle (Initial Puncture Needle - Three Part)

UROLOGY

INITIAL PUNCTURE NEEDLE - TWO PART

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking on the shaft

Advantage of Locking Hub IP Needle

During the procedure the locking hub helps in holding the inner needle with outer cannula together, which helps in easy initial puncture

Type

Two Part

Size (in G)

18

Length (in cm)

15 & 20



INITIAL PUNCTURE NEEDLE - TWO PART WITH ECHOGENIC TIP

Product Code	Size (in G)	Length (in cm)
BIPSO180152PB BIPSO180202PB	18 18	15 20

INITIAL PUNCTURE NEEDLE - TWO PART WITH ECHOGENIC TIP & LOCKING HUB

Product Code	Size (in G)	Length (in cm)
BIPSO18015LHB BIPSO18020LHB	18 18	15 20

SURE PUNCTURE™ NEEDLE

(INITIAL PUNCTURE NEEDLE - THREE PART)

Intended to make the initial percutaneous puncture to access the kidney and to enable the introduction of a guidewire during urological procedures (PCN, PCNL, Mini PERC & Micro PERC).

Key Features

- Siliconised needle for less trauma
- Echogenic tip for better visualization during ultrasound
- 1 cm marking

Advantage of Sure Puncture™ Needle

- Combination of 21 G Bevel tip and 16 G Cannula Needle
- Inner 21 G provides minimal site puncture and contrast injection to ensure desired puncture location
- Coaxial needle provides secured access
- Can introduce working guidewire with 21 G puncture

Type

Three Part

Size (in G)

21/16

Length (in cm)

20



SURE PUNCTURE™ NEEDLE WITH ECHOGENIC TIP

Product Code	Size (in G)	Length (in cm)
BSPSO160203PB	21/16	20

BIOPSY

- Biopsy Gun - Semi Automatic

UROLOGY

BIOPSY GUN - SEMI AUTOMATIC

Used for obtaining histological core samples from soft tissue and organs.

Key Features

- Echogenic and atraumatic needle tip
- Ergonomic and light weight
- Cannula with depth marking
- Easy and one hand operated with two stage arming
- Coaxial available with trocar and blunt needle

Biopsy Gun - Semi Automatic With Coaxial & Blunt Needle

- 14, 16 G x 10, 15 cm - For breast biopsy
- 16, 18 G x 15, 20 cm - For liver & kidney
- 18, 20 G x 15, 20 cm - For lung



BIOPSY GUN - SEMI AUTOMATIC WITH COAXIAL & BLUNT NEEDLE

Product Code	Size (in G)	Length (in cm)
BBG0014010BS2	14	10
BBG0014015BS2	14	15
BBG0016010BS2	16	10
BBG0016015BS2	16	15
BBG0016020BS2	16	20
BBG0018015BS2	18	15
BBG0018020BS2	18	20
BBG0020015BS2	20	15
BBG0020020BS2	20	20

ACCESSORIES

- Meatal Dilator - Male
- Dalela Female Self Urethral Dilator
- Bladder Evacuator
- Penile Clamp
- Irrigation Set (TUR Set)
- Cath Tap

UROLOGY

MEATAL DILATOR - MALE

Used for self-dilation of penile urethral meatus.

Key Features

- Designed specifically for dilation of meatal strictures
- Atraumatic design
- Uniform taper and smooth surface for easy insertion
- Superior handling grip

Size (in Fr)

6/30

MEATAL DILATOR - MALE

Product Code	Size (in Fr)
BMD0030000BM	6/30



DALELA FEMALE SELF URETHRAL DILATOR

PATENT PENDING

Used for self dilation of female urethra.

Key Features

- Novel design specifically for dilation of female urethra
- Made of durable, biocompatible material
- Soft tip for atraumatic and easy insertion
- Less trauma and pain during dilation
- Easy to use & extremely user friendly
- Uniform taper and smooth surface
- Can empty the bladder while device in place

Size (in Fr)

8/30

DALELA FEMALE SELF URETHRAL DILATOR

Product Code	Size (in Fr)
BMD0030000BD	8/30



BLADDER EVACUATOR

Used for irrigation and tissue collection during transurethral bladder and prostate surgery.

Key Features

- Used for TURP procedures
- Easy to fit connection to resectoscope instrument
- Soft easy to squeeze bottle design, allows fast and efficient removal of tissue during TURP

BLADDER EVACUATOR

Product Code

BBE0000000B1



PENILE CLAMP

Used for temporary control of male urinary incontinence.

PENILE CLAMP

Product Code

BPC0000000B0



IRRIGATION SET (TUR SET)

Used for endoscopic irrigation during transurethral resection of prostate gland.


Key Features

Non kinkable PVC tube

IRRIGATION SET - Y

Product Code

BTS0000000B0



CATH TAP

Used for intermittent drainage of foley catheter.

Key Features

- Intermittent drainage of foley catheter
- Reduce infection rate
- Avoid carrying of urine bags
- Can be fixed to patient thigh
- Useful for patients with long term catheterization

CATH TAP

Product Code

BCT0000000B0

